

	Stock solution			Working concentration	
	mg/ml:	solvent	store at	µg/ml	dilution
Ampicillin*	50	H <sub>2</sub> O	-20°C	20	0.4 µl/ml
Carbenicillin	50	H <sub>2</sub> O	-20°C	20	0.4 µl/ml
Chloramphenicol	34	EtOH	-20°C	25	0.7 µl/ml
Gentamycin	10	H <sub>2</sub> O	-20°C	15	1.5 µl/ml
Hygromycin B	10	H <sub>2</sub> O	-20°C	10	1 µl/ml
Kanamycin	50	H <sub>2</sub> O	+4°C	10	0.2 µl/ml
Kasugamycin	10	H <sub>2</sub> O	-20°C	1000	100 µl/ml
Neomycin	10	H <sub>2</sub> O	-20°C	100	10 µl/ml
Rifampicin*	34	EtOH	-20°C	150	4.4 µl/ml
Spectinomycin	10	H <sub>2</sub> O	-20°C	100	10 µl/ml
Streptomycin	50	H <sub>2</sub> O	-20°C	10	0.2 µl/ml
Tetracycline**	5-12.5	EtOH	-20°C	10	2 µl/ml

\* - Light sensitive;

\*\* - Light sensitive, Mg<sup>2+</sup> - inhibitor;